
Lehigh Integrated Nano-Characterization User Facilities

Linking to the University of Pennsylvania

Industrial Liaison Materials Research Cente Bethlehem, PA 18015-319 610-758-4855 http://www.lehigh.edu/~inmatso

Education Electronics Mechanical Structural Chemical



Nano-characterization in the Nano-Mechanical Behavior Laboratory

nning Probe Microscopy - Digital Instruments (DI)Nano surface imaging in air or in solution contact (AFM) or Tapping Mode electrochemical stages

AFM image of PTFE lines drawn out by

the friction transfer technique onto a glass

STM image of self-assembled Au

To analyze plastic deformation

materials at very low-load and

during shallow (only 3-30 nm) indentation

interpretation of nanoindentation

tests and to nano-asperity contact

Results are applicable to

Nano Indentation of a Platinum thin film with the

Analyze, model, and control stress relaxation

imensional stability for nm-scale metal films

measurement technique with unprecedented

by applying a new resonance based

Thin Film Creep Measurement by Membrane

Resonance. Operates up to 500 °C, with film thickness of

20-100 nm, stress resolution

typically < 1 MPa

Hysitron Triboscope

(150 nm "diameter", 14 nm depth)

nanoparticles on highly ordered

Contact and tapping modes, and electrochemical stages for control of reactions in -situ Imaging may be performed in air or under desired solution. A Hysitron Triboscope Nanoindentation attachment allows surface imaging and indentation, and measuremen of elastic and plastic properties, with indentation of 20 nm. Forces on order of uNto mN may be applied. Temperature control up to 150°C.

AFM Imaging of Nanostructures

Nano Indentation

Transmission Electron Microscopy

emission gun. It is equipped with energy dispersivespectroscopy (EDS) and an energy (EELS) and energy filtering of images and diffraction patterns.

of a bimodal array of self-organized gold in the IEOL 2010E

Analytical Electron Microscopy

Vacuum Generators HB 603 Scanning Transmission Electron Microscope.


This unique dedicated STEM operates up to 300kV. It is equipped with two EDS systems: one with a standard Si(Li) detector and one with an intrinsic Gedetector, Configured with an optimum geometry and large solid angles of collection. This makes our instrument the world leader for highresolution EDS and EELS analysis.

Lehigh Microscopy School

- > World renowned set of short courses in electron microscopy and related technique
- ➤ Over 5000 scientists, engineers and technicians from universities, industrial research labs and national labs trained
- > U.S. experts educate participants in scanning, transmission, analytical and othermicroscopies and spectroscopies

- Atomic Force and Scanning Probe
- Microscopy (AFM/SPM)
- *Analytical Electron Microscopy (AEM)

Lattice Imaging of Nanoparticles



(b) simulated HREM image (c) structural model

A unique high resolution, high

sensitivity instrument for analysis

of the nano chemistry of surfaces The Scienta has a rotating anode and monochromatized electron gun with a low energy electron diffraction/residual gas analyzer. It operates in the range -175 C to >

600 °C, UHV to one atmosphe

AEM of Nanostructures

This X -ray map shows Yttrium segregation to grain ndaries in a Y -doped ZrO2 ceran The field of view is 300 x 300 nm2





STEM-XEDS elemental maps from Au-Ag core-shell nanoparticles

Surface Chemical Analysis- Angle Resolved XPS

Scienta ESCA-300







Nanoelectronics Educational Programs





EXTENSIVE FACILITIES









Internet2 videoconferencing for shared courses and research collaboration

Managing the Impact of Nanotechnology of

Society and the Environment

 Shared research facilities and instrumentation between PA universities Smartboard & Web-based course materials for shared exchange of

chigh maintains the nation's 3rd oldest Science, Technology, & Society Program year tradition of collaborative interdisciplinary research in the social, ethical, an political implications of scientific research and technological innovation

Science, Nanotechnology, & Society Programs are needed :

- To sensitize people as early as possible to the anticipated, and profound social consequences of nanotechnology, especially in the area of molecu lar biotechnology, defense, and robotics > To encourage cross-disciplinary thinking to accomplish constructive, rather that
- destructive, nanotechnology applications > To consider the social and environmental consequences of nanotechnology as r
- > To promote the positive impacts of nanotechnology on the environment (eg improved sensing techniques and processes for removal of the sma llest contaminants in water soil and air)
- > To debunk misinformation and disseminate accurate information or

Thin Film Creep Studies Optical Spectroscopy and Luminescence Microscopy of Nano-Structures

Focused Ion Beam Thinning

FEI Strata DB 235

This Dual-Beam FIB is a combination of a

focused ion beam system (FIB). It is used for

preparing TEM samples from specific localized

We combine the spatial resolution of a confocal microscope and n ear field optical microscope with the high spectral resolution and selectivity offered by combined excitation - emission spectroscopy. This allows correlation of atomic scale defects with optical properties observed on thenanoscale

Stages of preparation of piezoelectric ceramic TEM-foil



Domains in a confocal luminescence microscope (contrast is obtained through changes in the atomistic defect structure which also determine mobility, shape of the domains and the voltage required for domain switching)

Optical Molecular Spectroscopies

Infra Red Spectroscopy (Complementary to Raman) UV-VIS Spectroscopy (Oxidation States of Metal Oxides, as small as 0.4 nm)

In-situ Molecular Characterization of Nano-

Oxides Under Different Temperature.

Pressure and Gaseous/Liquid Environments

Laser Raman Spectroscopy (Structure of Metal Oxides, as small as 0.4 nm)

T = -200 to 1500°C